

Running the eSPC tools on Windows

February 2026



1. Install WSL2

- a. Open **PowerShell as Administrator** and run:
\$ wsl --install
- b. Verify the installation and restart the PC if required:
\$ wsl --status

2. Install a Linux distribution

- a. Open **PowerShell as Administrator** and run:
\$ wsl --install -d Ubuntu

3. Install Docker

- a. Go to the official Docker website and follow the instructions to install the Docker engine:
<https://docs.docker.com/engine/install/ubuntu/>

4. Pull the docker image and run it

- a. In the Ubuntu **Terminal** run:
 - i. For **ThermoAffinity**:
\$ docker pull emblspc/thermoaffinity_espc:1.0
\$ docker run -p 3838:3838 emblspc/thermoaffinity_espc:1.0
 - ii. For **MoltenProt**:
\$ docker pull emblspc/moltenprot_espc:1.1
\$ docker run -p 3838:3838 emblspc/moltenprot_espc:1.1
 - iii. For **FoldAffinity**:
\$ docker pull emblspc/foldaffinity_espc:1.0
\$ docker run -p 3838:3838 emblspc/foldaffinity_espc:1.0
 - iv. For **ChiraKit**:
\$ docker pull emblspc/chirakit_espc:1.0
\$ docker run -p 3838:3838 emblspc/chirakit_espc:1.0
 - v. For **PhotoMol**:
\$ docker pull emblspc/photomol_espc:2.0
\$ docker run -p 3838:3838 emblspc/photomol_espc:2.0
 - vi. For **KinGenie**:

```
$ docker pull emblspc/kingenie_espc:1.0  
$ docker run -p 3838:3838 emblspc/kingenie_espc:1.0
```


vii. For **Raynals**:

```
$ docker pull emblspc/raynals_espc:1.0  
$ docker run -p 3838:3838 emblspc/raynals_espc:1.0
```

5. Open the tool

- a. Open any Browser and go to <http://0.0.0.0:3838/>

Congratulations! The tool should be open and ready for use.



Condition	Include	Read	Σ
carbonic anhydrase A1 12:02:17	<input checked="" type="checkbox"/>	1	
carbonic anhydrase A2 12:05:51	<input checked="" type="checkbox"/>	1	